

GENTHERM

Astoflo® Plus Eco Blood and Infusion Warmer



Astoflo Plus Eco is a universal blood and infusion warmer that does not require any special disposables. Fluids are warmed by placing the standard infusion set in the heating profile and are kept warm all the way to the patient cannula.

Features and Benefits




The electric heating profile with integrated temperature sensors guarantees optimum warming and maximum safety.



The heating profile actively warms the fluid as it passes along the infusion line, and will also maintain the temperature of any prewarmed solutions, right to the patient's cannula.



The Astoflo provides an easy and economical alternative to warming intravenous infusions.



This consumable-free device helps reduce single-use plastic in healthcare settings.

Astoflo® Ordering Information

Product Description	1pcs / box Product Code:
ASTOFLO PLUS ECO	
Astoflo Plus Eco Blood and Infusion Warmer	AFP300
Heating profile for Astoflo Plus Eco warmer 4-5mm OD x 180cm length	WP31

