

GENTHERM

Electri-Cool® II Localized Cold Therapy



HELP REDUCE SWELLING AND PAIN

The Electri-Cool II delivers convenient and fast cooling therapy to help reduce inflammation and pain.

This portable, easy-to-use and lightweight system requires no ice to operate; simply add sterile water, connect a Cold Therapy Pad and select your preferred temperature to initiate therapy.



Portable, lightweight and user-friendly.



Three temperature settings designed for personalized comfort.



Helps to reduce swelling and associated pain.

Perfect for Orthopedics, Plastics, Dental and Oral Surgery, General Surgery & more! Enhance post-op recovery and provide more pain management options with Electri-Cool II.

No ice required!

Avoid the hassle, limitations and potential dangers of using traditional cooling methods such as ice:

- Potential skin damage from prolonged or direct application of ice to the skin.¹
- Condensation and rapid temperature loss.
- No temperature control.

Consistent, Active Cooling with Cold Therapy Pads & Wraps

Gentherm Cold Therapy Pads are designed to deliver precise, controlled cold therapy. Paired with the Electri-Cool II, the Pads maintain a consistent temperature and remain pliable to conform to the contours of the body, even at the coldest temperatures.



1. Cleveland Clinic medical. (2023, February 22). Frostbite: Signs & symptoms, stages, treatment & prevention. Cleveland Clinic. <https://my.clevelandclinic.org/health/diseases/15439-frostbite>

PRODUCT SPECIFICATIONS

Dimensions:

12" W x 6" D x 12" H
(30.5 cm x 15.2 cm x 30.5 cm)

Weight:

13 lbs. (5.9 kg)

Water Temperature Range:

5°C-7°C (40°F-45°F)
8°C-10°C (45°F-50°F)
11°C-13°C (50°F-55°F)

Electri-Cool® II System Ordering Information

CAT #	Product Description	Voltage/Quantity
ELECTRI-COOL II		
767	Electri-Cool II	100v / 115v, 220v / 240v
757 H1D	12" (30.5 cm) Dual Pad Connecting Hose	Ea.
757 HIE	6' Connecting Hose (Insulated for Electri-Cool II Model 757)	Ea.
767 -BBK	Bed Bracket	Ea.
UNV LPS	Universal Stand	Ea.
EC -080	Adapter Plate (for use with universal stand)	Ea.
TEMP-PADS® COLD THERAPY PADS		
814-NS	18.5" x 13" (47cm x 33cm) with straps, Non-Sterile	10 / box
1112-NS	11" x 12" (27.9cm x 30.5cm) with straps, Non-Sterile	10 / box
812-NS	8" x 12" (20.3cm 30.5cm) with straps, Non-Sterile	10 / box
CT-99 NS	10" x 11" (25.4cm x 28cm) Shoulder/Knee with straps, Non-Sterile	10 / box




CT-99 NS



1112-NS



814-NS



812-NS