

# GENTHERM

## Micro-Temp® LT Localized Heat Therapy



### REDUCE PAIN AND PROMOTE HEALING

Help reduce pain, promote comfort and kick-start rehabilitation with the Micro-Temp LT. User-friendly and portable, Micro-Temp LT is designed to deliver targeted warmth directly to areas of discomfort or pain.

Set the water temperature between 20°C and 42°C (68°F - 107.6°F), apply a specialized Gentherm Maxi-Therm® Lite pad, and be confident your patient is receiving the optimal benefits of heat therapy.



Improve blood flow to injured areas



Reduce pain from muscle spasms and stiffness



Relief from discomfort



Help with the rehabilitation from medical procedures

*Perfect for the Emergency Department, MedSurg, Labor & Delivery, Mother/Infant and more!*

## PRODUCT SPECIFICATIONS

### Dimensions:

9.8" W x 8.2" D x 9.1" H  
(24.9 cm x 20.8 cm x 23.1 cm)

### Weight:


9 lbs (4.1 kg)

### Water Temperature Range:

20°C-42°C (68°F-107.6°F)

# Micro-Temp® LT System Ordering Information

CAT #	Product Description	Voltage/Quantity
<b>MICRO-TEMP LT</b>		
749	Micro-Temp LT	100v / 115v, 220v / 240v
UNV LPS	Universal Stand	Ea.
MT-081	Adapter Plate (use with universal stand)	Ea.
<b>MAXI-THERM® LITE PADS (SINGLE-PATIENT USE)</b>		
871	25" x 17.25" (63.5cm x 43.82cm)	10 / box
872	25" x 19" (63.5cm x 48.26cm)	10 / box



Maxi-Therm Lite Pads