

DAYMARK™ | SURGICAL

PORTABLE LED HEADLIGHT SYSTEM



SEE BETTER
SEE FARTHER
SEE WHY




**SINCE
1971** CELEBRATING
FIVE DECADES
OF ILLUMINATION

DAYMARK™ SURGICAL

HIGH-INTENSITY, PORTABLE LED SURGICAL HEADLIGHT WITH ADJUSTABLE SPOT, POWERED BY A 4-CELL BATTERY

- Illumination Range: ~60,000 LUX at 16" (40 cm) to ~260,000 LUX at 10" (25 cm)
- Adjustable Spot Diameter: 1"-4" (2.5 - 10.16 cm) at a working distance of 16" (40 cm)
- Color Temperature: 5,000 K
- Intensity Settings: 5 Levels
- 4-Cell Battery (System includes 2 batteries)
- Battery Duration at Full Intensity > 5 hours
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 50,000 hours
- Optic Weight: 5 oz (141 g)
- Battery Weight: 8.5 oz (240 g)



DAYMARK PORTABLE LED HEADLIGHT SYSTEM INCLUDES:

- Daymark Surgical Headlight on Fiji Headband
- Two 4-Cell Batteries
- Charging Cradle, Power Supply
- Holster, Linkage Kit

DAYMARK PORTABLE LED HEADLIGHT SYSTEM PRODUCT CODES:

- | | |
|-------------|---|
| BFW 9070 | Daymark Headlight System with Fiji Headband with Two 4-Cell Batteries |
| BFW 5440 | Daymark Protective Carrying Case, molded interior provides space for all components |
| BFW 9022ACC | Daymark Mini Battery Charger |
| BFW 9015 | Daymark 4-Cell Battery |
| BFW 9076 | Daymark Holster |

EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS

- | | | | | |
|--------------------|--------------------|-----------------------|--------------------|--------------------|
| HOSPITAL | ENT | SURGERY CENTER | ORTHOPEDIC SURGERY | |
| | GENERAL SURGERY | | PLASTIC SURGERY | |
| | CARDIOSURGERY | | | |
| | NEUROSURGERY | | OFFICE | DENTAL EXAM |
| | ORTHOPEDIC SURGERY | | | ORAL/MAXILLOFACIAL |
| PEDIATRIC SURGERY | | | | |
| SPINE SURGERY | | | | |
| THORACIC SURGERY | | | | |
| TRANSPLANT SURGERY | | | | |
| VASCULAR SURGERY | | | | |



BFWinc.com

800-717-4673 | +1-502-899-1808

1550 SAMCO ROAD
RAPID CITY, SD, 57702 USA

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2024 BFW, Inc.

ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE

DAYMARK TESTED AND CERTIFIED TO: IEC 60601-1 Edition 3, IEC 60601-1-2 Edition 4, ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10), CAN/CSA C22.2 NO 60601-1 Edition 2 (2008), CE